## Amber<sup>®</sup> Mill Instructions for Use

#### 1. Overview

- · Trade / Device Name : Amber Mill
- · CommonName : Dental Frame Material for Dental Prosthesis
- Intended Use of the Device : Amber Mill Series are indicated for fabricating glass ceramic restorations such as single-unit anterior and posterior crowns, veneers, inlays/onlays, and anterior 3-unit bridges using CAD/CAM system.
- · Classification Name : Porcelain Powder for Clinical Use
- · Packaging Unit : Refer to HASS standard package.

#### 2. Instruction



(1) How to Use

This product must be used in accordance with the using methods of a dental CAD/CAM system.

\* Procedure for using glass ceramic blocks \*

- ① Attach the Jig to the accurate location.
- Mount it on the CAD/CAM equipment.
- ③ Inputs the size information of the prepared block into the CAD/ CAM equipment.
- ④ Inputs correction information needed for processing
- ⑤ Process the Block using the CAD/CAM equipment.
- 6 Carefully detach the process-completed block from the equipment.
- ⑦ Detach the processed artificial tooth or restoration from the block.
  ⑧ Artificial teeth or prosthetic separated by heat treatment at 810 ~
- 865°C makes
- crystallization.
- (9) If necessary, perform stain and glazing treatment.
- (2) Storage and Maintenance before Use

Do not store in package open or dirty place it may contaminate the products.

- ② Store away from moisture, direct sunlight, and heat.
- ③ Do not reuse or recycle the remaining part once used.

#### $\triangle$ 3. Cautions

- (1) Cautions before Use
  - Be careful not to damage the milling tool of the CAD/CAM machine when attaching
    - or detaching the product.
  - ② Be careful not to get your hand caught in the milling tool.
  - ③ The jig should be attaching to an accurate location.
  - ④ Suppress or remove the dust which may occur during the operation of CAD/CAM machine.
  - $\ensuremath{\textcircled{}}$  Do not drop the product on the ground or apply heavy force as it may damage the
    - product.
  - 6 Keep the product out of reach of infants and children.
  - ⑦ Product should be handled by dental technicians and dentists.
- (2) Storage and Maintenance before Use
  - Store the product at room temperature in a dry place.
  - O Pack and store the product properly to ensure that it is not damaged.
  - ③ Store the product at temperatures ranging from 0°C ~ 40°C, in combination with relative humidity of 10% r.H ~ 90% r.H, under atmospheric pressures ranging from 500 hPa ~ 1060 hPa.

#### 4. Side effect

It the patient is known to be allergic to any of the components of Amber Mill, the material must not be used to fabricate restorations.

#### 5. Contraindication

- ① Posterior bridges reaching into the molar region
- ② 4-and more-unit bridges
- ③ Inlay-retained bridges
- ④ Very deep sub gingival preparations
- ⑤ Bruxism
- 6 Cantilever bridges / extension units
- ⑦ Maryland bridges
- ③ Any other use not listed in the indications

#### 6. Mechanical and Physical Properties

- ① Material : Glass-ceramics
- ② Flexural Strength : over 300 MPa
- ③ Chemical Solubility : below 100 µg/cm²
- ④ Coefficient of Thermal Expansion : 10.0 (±0.5) x 10<sup>-6</sup> K<sup>-1</sup>
- \* This is a single-use product. \* Do not reuse.

#### 7. Pictograph



#### <€2195

KTR Europe GmbH Mergenthalerallee 77, 65760 Eschborn, Germany

#### HASS Corporation

77-14, Gwahakdanji-ro, Gangneung-si, Gangwon-do, KOREA 25452 Tel: +82-70-7712-1300 / Fax: +82-33-644-1231 Customer Support : +82-2-2083-1367 E-mail : hasscorp@hassbio.com Website : www.hassbio.com



# Amber<sup>®</sup> Mill Recommended Toughening Heat-treatment Schedule

It is possible to differentiate translucency with a single block of Amber<sup>®</sup> Mill. Just decide what shade you will use, then choose the toughening heat-treatment temperature according to your targeted translucency. This will enhance the efficiency in work process and inventory management for CAD/CAM milling blocks.

## VITA VACUMAT

Predry ℃	 min.	 min.			⊤ ℃		→ min.	VAC min.		°*
400	3.00	HT	6.50	60	ΗT	815	15.00	ΗT	21.50	690
		MT	7.05		MT	825		MT	22.05	
		LT	7.20		LT	840		LT	22.20	
		MO	7.40		MO	860		MO	22.40	

\* The firing chamber must not be opened during long term cooling.

VACUMAT is a registered trademark of VITA.

## IVOCLAR VIVADENT PROGRAMAT CS

B °C	S min.	t ≯ ℃ / min.	T °C	;	H min.	H VAC. 1 °C / in. VAC. 2 °C		L C	tL∗	
400 3.		60	HT	815	15.00	HT	550/815	690	0	
	3.00		MT	825		MT	550/825			
			LT	840		LT	550/840			
			MO	860		MO	550/860			

\* The firing chamber must not be opened during long term cooling.

PROGRAMAT CS is a registered trademark of IVOCLAR VIVADENT.

## DEKEMA Austromat 624i

	HT				MT			LT(31:10)			MO		
Dry			-;			-:-						-;	
Close			02:00			02:00			02:00			02:00	
Preheat	450°C		01:00										
Temperature 1	830°C	60°C /min	15:00	840°C	60°C /min	15:00	855°C	60°C /min	15:00	875°C	60°C /min	15:00	
Temperature 2	690°C	60°C /min	-;	690°C	60°C /min	-;	690°C	60°C /min	-÷	690°C	60°C /min	;	
Temperature 3	°C	°C/ min	-;	°C	°C/ min	-:-	°C	°C/ min	;	°C	°C/ min	;	
VAC (off/level/hold)	830°C	100%	15:00	840°C	100%	15:00	855°C	100%	15:00	875°C	100%	15:00	

\* The firing chamber must not be opened during long term cooling.

Austromat 624i is a registered trademark of DEKEMA.

