# HeraCeram Zirkonia 750

|      |                                  |    |
|------|----------------------------------|----|
| (DE) | Gebrauchsanweisung .....         | 3  |
| (GB) | Instructions for use .....       | 6  |
| (FR) | Mode d'emploi .....              | 9  |
| (ES) | Instrucciones de uso .....       | 12 |
| (IT) | Istruzioni per l'uso .....       | 15 |
| (PT) | Instruções de uso .....          | 18 |
| (NL) | Gebruiksaanwijzing .....         | 21 |
| (SE) | Bruksanvisning .....             | 24 |
| (DK) | Brugervejledning .....           | 27 |
| (NO) | Bruksanvisning .....             | 30 |
| (FI) | Käyttöohjeet .....               | 33 |
| (GR) | Οδηγίες χρήσης .....             | 36 |
| (CZ) | Návod k použití .....            | 39 |
| (HU) | Használati utasítás .....        | 42 |
| (LV) | Lietošanas instrukcija .....     | 45 |
| (LT) | Naudojimo instrukcija .....      | 48 |
| (PL) | Instrukcja obsługi .....         | 51 |
| (HR) | Upute za uporabu .....           | 54 |
| (RO) | Instrucţiuni de utilizare .....  | 57 |
| (TR) | Kullanma talimatı .....          | 60 |
| (RU) | Инструкция по применению .....   | 63 |
| (AE) | تعليمات الاستخدام .....          | 66 |
| (CN) | 使用说明 .....                       | 69 |
| (KR) | 사용 설명서 .....                     | 72 |
| (SK) | Návod na použitie .....          | 75 |
| (BG) | Инструкции за употреба .....     | 78 |
| (UA) | Інструкція із застосування ..... | 81 |
| (SI) | Navodila za uporabo .....        | 84 |
| (RS) | Uputstva za upotrebu .....       | 87 |
| (MK) | Упатство за употреба .....       | 90 |

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Stand: 2024-03

**Medical Device – for use by dental health care professionals only.**

## Intended purpose

Dental ceramic for veneering, customisation and modification of prosthetic restorations.

Classification according to EN ISO 6872: Type I; Class 1

This dental material is suitable for patients requiring dental treatment for the following indications with consideration of the contraindications. For pregnant and nursing women, the treatment risks must be weighed carefully against the benefits taking into consideration the unborn child or infant.

Kulzer medical devices ensure the rehabilitation of oral functions such as chewing, speaking and aesthetics. They stabilize the remaining dentition and/or the alveolar ridge restoratively or prosthetically.

## Indications

Veneering material for ceramic frameworks, e.g. manufactured from zirconium dioxide and lithium disilicate, with a CTE/ coefficient of thermal expansion ( $\alpha_{25-500^{\circ}\text{C}}$ ) of 9.7 – 11.0  $\mu\text{m/m-K}$ .

## Contraindications

- Bruxism or other parafunctional habits.
- The use of this product is contraindicated in case of known or suspected allergies against components of this product.

## Potential side effects

Hypersensitivities to the product or its components cannot be excluded in individual cases.

## Composition

### HeraCeram Zirkonia 750 powder materials

Glass ceramic, chromium-tin-pink-sphene, vanadium-zirconium-yellow-baddeleyite, diyttrium oxide silicate, tin(IV) oxide

### HeraCeram Zirkonia 750 Adhesive (paste)

Glass ceramic, 1,4-butanediol, urea

## Application notes

- HeraCeram Zirkonia 750 may only be used on indicated framework materials whose coefficient of thermal expansion (CTE) is within the approved CTE range.
- HeraCeram Zirkonia 750 must not be mixed/combined with veneering ceramic components of ceramic lines manufactured by competitors.
- The use of modelling/staining liquids or thinner liquids of ceramic lines manufactured by competitors is not indicated.
- Powdery materials (e.g. dentine, incisal, etc.) must not be mixed with paste-like materials (e.g. adhesive, stains, glaze).
- If the consistency of the HeraCeram Zirkonia 750 Adhesive paste becomes drier and thus firmer due to prolonged storage, the viscosity can be readjusted by carefully adding the HeraCeram Paste Opaque Liquid (POL).
- Frameworks for ceramic veneering should have the reduced anatomical shape. Sharp edges and corners on the frameworks should be avoided.

## Product application

### Veneering of zirconium dioxide frameworks

#### 1. Framework preparation

After completion of framework processing, the items must be thoroughly cleaned of abrasive dust and abrasive residues. The zirconium oxide frameworks must be cleaned with the steam jet. It is not necessary to blast the veneer surface to increase the adhesive bond.

#### 2. Adhesive

For an optimum bond between zirconium dioxide and HeraCeram Zirkonia 750, HeraCeram Zirkonia 750 Adhesive is applied to the framework surface with a brush in an evenly thin layer (wash) and fired under vacuum (see firing chart).

#### 3. Layering

The ceramic veneer is designed in a two-layer technique (dentine incisal layering) or, depending on the desired aesthetic requirements, with additional individual materials<sup>2</sup>, and is fired. The margin ceramics LM and correction material (COR) are used for corrections after completing the veneering. The margin ceramics HM are used to produce all-ceramic crown margins. To produce the ready-to-use state, the ceramic powders are mixed with the modelling liquids ML or MLS from Kulzer's HeraCeram range. HeraCeram Zirkonia 750 margin ceramics HM and LM are mixed with HeraCeram margin ceramic liquid (SML). For the firing procedure, see the firing table.

**Note:** To prevent breathing in the ceramic dust when grinding the ceramic materials, wear a face mask and eye protection and work with a suction device.

#### 4. Glaze firing

The veneer can be individually characterised with HeraCeram Glaze universal and HeraCeram Stains universal stains according to the Instructions for Use. The degree of glazing and the texture of the ceramic surface are defined during glaze firing by the temperature and the holding time until the final temperature. For this reason, the specifications for the glaze firing temperature and holding time are for orientation only and may need to be adjusted to achieve the required result. Firing procedure see firing table. It is also possible to mechanically polish HeraCeram Zirkonia 750.

### Veneering of lithium di-silicate frameworks

#### 1. Framework preparation

After completion of framework processing, the items must be thoroughly cleaned of abrasive dust and abrasive residues, e.g. using a steam jet. Do not blast the frame with  $\text{Al}_2\text{O}_3$  or glass beads!

#### 2. Adhesive

For an optimum bond between lithium di-silicate and HeraCeram Zirkonia 750, Adhesive 750 is subsequently applied to the framework surface with a brush in an evenly thin layer (wash) and fired under vacuum (see firing chart).

#### 3. Layering

See above. Please pay attention to the firing procedure.

#### 4. Glaze firing

See above. Please pay attention to the firing procedure.

Colour mapping table

| Shade                                | A1    | A2    | A3    | A3,5    | A4    | B1    | B2    | B3    | B4    | C1    | C2    | C3    | C4    | D2    | D3    | D4    | BL1   | BL3   |
|--------------------------------------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| Increaser                            | IN A1 | IN A2 | IN A3 | IN A3,5 | IN A4 | IN B1 | IN B2 | IN B3 | IN B4 | IN C1 | IN C2 | IN C3 | IN C4 | IN D2 | IN D3 | IN D4 |       |       |
| Dentine                              | D A1  | D A2  | D A3  | D A3,5  | D A4  | D B1  | D B2  | D B3  | D B4  | D C1  | D C2  | D C3  | D C4  | D D2  | D D3  | D D4  | D BL1 | D BL3 |
| Chromadentine                        | CD A1 | CD A2 | CD A3 | CD A3,5 | CD A4 | CD B1 | CD B2 | CD B3 | CD B4 | CD C1 | CD C2 | CD C3 | CD C4 | CD D2 | CD D3 | CD D4 |       |       |
| Incisal                              | S1    | S1    | S2    | S2      | S4    | S1    | S1    | S2    | S4    | S1    | S3    | S3    | S3    | S1    | S2    | S2    | S BL  | S BL  |
| Mamelon-dentine<br>Secondary dentine | MD1   | MD1   | SD2   | SD2     | SD2   | MD2   | MD2   | MD3   | MD3   | MD2   | SD1   | SD2   | SD2   | MD1   | MD3   | SD1   |       |       |
| Value                                | VL1   | VL2   | VL3   | VL4     | VL4   | VL1   | VL2   | VL3   | VL4   | VL1   | VL2   | VL3   | VL4   | VL2   | VL3   | VL4   |       |       |
| Opal incisal                         | OS1   | OS1   | OS2   | OS2     | OS4   | OS1   | OS1   | OS2   | OS4   | OS1   | OS3   | OS3   | OS3   | OS1   | OS2   | OS2   |       |       |
| Margin ceramics<br>HM/LM             | 1     | 1     | 2     | 2       | 6     | 3     | 3     | 4     | 4     | 1     | 1     | 6     | 6     | 1     | 2     | 4     |       |       |

General firing table

**IMPORTANT NOTICE:** The firing temperatures given are intended as guidelines. Deviations are possible due to varying furnace power and may need to be adapted.

| General firing table<br>zirconium dioxide |        | Adhesive<br>750 | Margin<br>ceramics HM | 1. Dentine <sup>2</sup> | 2. Dentine <sup>2</sup> | Stains-/ Glaze<br>universal | Correction<br>material COR | Margin<br>ceramics LM |
|-------------------------------------------|--------|-----------------|-----------------------|-------------------------|-------------------------|-----------------------------|----------------------------|-----------------------|
| Start temperature                         | °C     | 500             | 500                   | 500                     | 500                     | 500 <sup>1</sup>            | 500                        | 500                   |
| Drying                                    | min    | 5               | 5                     | 5                       | 5                       | 5                           | 4                          | 4                     |
| Pre-heating                               | min    | 1               | 2                     | 2                       | 2                       | 1                           | 2                          | 2                     |
| Temperature increase                      | °C/min | 75              | 50                    | 50                      | 50                      | 60                          | 50                         | 50                    |
| Final temperature                         | °C     | 800             | 780                   | 750                     | 745                     | 725                         | 715                        | 715                   |
| Holding time                              | min    | 2               | 1                     | 1                       | 1                       | 0,5-1                       | 1                          | 1                     |
| Vacuum start                              | °C     | 500             | 500                   | 500                     | 500                     | –                           | 500                        | 500                   |
| Vacuum stop                               | °C     | 800             | 780                   | 750                     | 745                     | –                           | 715                        | 715                   |

| General firing table<br>lithium di-silicate |        | Adhesive<br>750 | Margin<br>ceramics HM | 1. Dentine <sup>2</sup> | 2. Dentine <sup>2</sup> | Stains-/ Glaze<br>universal | Correction<br>material COR | Margin<br>ceramics LM |
|---------------------------------------------|--------|-----------------|-----------------------|-------------------------|-------------------------|-----------------------------|----------------------------|-----------------------|
| Start temperature                           | °C     | 400             | 400                   | 400                     | 400                     | 400 <sup>1</sup>            | 400                        | 400                   |
| Drying                                      | min    | 5               | 5                     | 5                       | 5                       | 5                           | 4                          | 4                     |
| Pre-heating                                 | min    | 1               | 2                     | 2                       | 2                       | 1                           | 2                          | 2                     |
| Temperature increase                        | °C/min | 50              | 50                    | 50                      | 50                      | 60                          | 50                         | 50                    |
| Final temperature                           | °C     | 750             | 780                   | 750                     | 745                     | 725                         | 715                        | 715                   |
| Holding time                                | min    | 1               | 1                     | 1                       | 1                       | 0,5-1                       | 1                          | 1                     |
| Vacuum start                                | °C     | 400             | 400                   | 400                     | 400                     | –                           | 400                        | 400                   |
| Vacuum stop                                 | °C     | 750             | 780                   | 750                     | 745                     | –                           | 715                        | 715                   |

Cooling rate: There are no defined cooling rates required for HeraCeram Zirkonia 750 veneering ceramics.

<sup>1</sup> Start temperature 600°C for stain shade and glaze pastes!

<sup>2</sup> implies the following masses: Dentine / Chroma dentine / Incisal / Transpa / Mamelondentine / Secondary dentine / Value / Opal incisal / Opal transpa / Enhancer / Mask / Gingiva / Increaser

**Technical data**

**Dentine / Chroma dentine / Incisal / Transpa / Mamelondentine / Secondary dentine / Value / Opal incisal / Opal transpa / Enhancer / Mask / Gingiva / Increaser:**

|                |                  |
|----------------|------------------|
| CTE (25-400°C) | 9.0 ± 0.2 µm/m-K |
| CTE (25-500°C) | not available    |
| TG             | 480 ± 20 °C      |

**Correction material COR / Margin ceramics LM:**

|                |                  |
|----------------|------------------|
| CTE (25-400°C) | 8.9 ± 0.2 µm/m-K |
| CTE (25-500°C) | not available    |
| TG             | 470 ± 20 °C      |

**Margin ceramics HM:**

|                |                  |
|----------------|------------------|
| CTE (25-400°C) | 9.0 ± 0.2 µm/m-K |
| CTE (25-500°C) | not available    |
| TG             | 510 ± 20 °C      |

**Adhesive 750:**

|                |                  |
|----------------|------------------|
| CTE (25-400°C) | 9.0 ± 0.2 µm/m-K |
| CTE (25-500°C) | not available    |
| TG             | 485 ± 20 °C      |

If the glass transition temperature is approximately 500 ± 30°C and there is a significant change in slope of the CTE curve before and after the glass transition temperature, then the value of 400°C is stated.

**General information**

Please pay attention to the following safety instructions and all processing instructions in other sections. This product must be used according to the instruction for use only. Any use not complying with these instructions is at the discretion and sole responsibility of the dental professional itself. Visually inspect for damage before use. Damaged products must not be used.

**Warnings and safety instructions**

Avoid breathing dust/fume/gas/mist/vapours/spray. Ensure adequate ventilation and exhaust. Wear suitable face mask and eye protection. Inhalation can cause health problems. May cause drowsiness or dizziness. These dangers and safety instructions do not refer to the prosthetic restoration itself, but to the vapours and dusts released during the processing of the material.

**Special advice**

For SSCP please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact [sscp-service@kulzer-dental.com](mailto:sscp-service@kulzer-dental.com)  
To find the SSCP in the EUDAMED use the UDI code attached at the packaging of the product.

**Storage conditions**

Carefully reseal the container after removing the product to avoid contamination, clumping and drying out. No special storage conditions required.

**Disposal information**

Recommendation: Dispose of in accordance with official regulations. Do not dispose of contents or partially emptied packaging in the household waste or allow them to enter the sewage system.  
European Waste Catalogue: 180106 Chemicals consisting of or containing hazardous substances.

**Complaint notice**

Please quote batch number [LOT] and article number [REF] or UDI in all correspondence about the product.  
According to EU Medical Device Regulation users/ patients are obliged to report serious incidents with a medical device to the manufacturer and to the competent authority of the country, the incident occurred.

Safety data sheets and more information are available at our website [www.kulzer.com](http://www.kulzer.com)

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Kulzer GmbH